

Anti-Myh4/Myh2B Rabbit Monoclonal Antibody For Wb - Q9Y623

Item Number: CM0011555

Introduction

High-quality rabbit monoclonal antibody targeting myosin-4 (MYH4/MyHC-IIb) for WB, IF-P, IHC-P, ELISA. Cross-reacts with human, mouse, rat. 223kDa. Ideal for muscle research.

[Learn More](#)

Parameter	Value
Product Name	MYH4 Rabbit mAb
Remarks/Alias	MYH2B; MyHC-2B; MyHC-IIb
Species	Human
GeneID (Human)	4622
GeneID	4622
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1500-1600 of human MYH4 (NP_060003.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y623
Protein Weight	223kDa
Shipping	Ice bag