


Anti-Myh6/A-Mhc Rabbit Polyclonal Antibody - Q02566

Item Number: CM0026875



Introduction

Rabbit polyclonal antibody against mouse MYH6/α-MHC (Myosin-6) for WB, IP, IF/ICC, IHC-P, ELISA. Cross-reacts with mouse and rat. Ideal for cardiac muscle research.

[Learn More](#)

Parameter	Value
Product Name	MYH6/α-MHC Rabbit pAb
Remarks/Alias	Myhca; Myhc-a; alphaMHC; alpha-MHC; A830009F23Rik
Species	Mouse
GeneID (Human)	4624
GeneID	17888
Immunogen	Synthetic peptide
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IP, IF/ICC, IHC-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q02566
Protein Weight	224kDa
Shipping	Ice bag