

## Anti-Myh6/Myh7 Monoclonal Antibody For Wb, If-P, Elisa - P13533 / P12883

Item Number: CM0008754

### Introduction

MYH6/MYH7 Rabbit mAb is a monoclonal antibody validated for Western blot, immunofluorescence, and ELISA applications. Detects human, mouse, and rat cardiac myosin heavy chain isoforms. Ideal for cardiovascular disease research.

[Learn More](#)

| Parameter        | Description                                                                                                                      |
|------------------|----------------------------------------------------------------------------------------------------------------------------------|
| Product Name     | MYH6/MYH7 Rabbit mAb                                                                                                             |
| Remarks/Alias    | ASD3; MYHC; SSS3; CMH14; MYHCA; CMD1EE; alpha-MHC; CMH1; MPD1; SPMD; SPMM; CMD1S; MYHCB                                          |
| Species          | Human                                                                                                                            |
| Gene ID (Human)  | 4624/4625                                                                                                                        |
| Immunogen        | Synthetic peptide: A synthetic peptide corresponding to a sequence within amino acids 1730-1830 of human MYH6/MYH7 (NP_002462.2) |
| Source           | Rabbit                                                                                                                           |
| Category         | Monoclonal Antibodies                                                                                                            |
| Application      | WB, IF-P, ELISA                                                                                                                  |
| Cross Reactivity | Human, Mouse, Rat                                                                                                                |
| SWISS            | P13533/P12883                                                                                                                    |
| Protein Weight   | 223kDa                                                                                                                           |
| Shipping         | Ice bag                                                                                                                          |