

Anti-Myh7 Rabbit Polyclonal Antibody - P12883

Item Number: CM0024256



Introduction

Rabbit polyclonal antibody targeting human MYH7 (myosin-7), a cardiac/slow skeletal muscle motor protein. Validated for WB, IF-P, IHC-P, ELISA; cross-reacts with human, mouse, rat. Suitable for muscle contraction and cardiomyopathy studies.

[Learn More](#)

Parameter	Value
Product Name	MYH7 Rabbit pAb
Remarks/Alias	CMH1, MPD1, SPMD, SPMM, CMD1S, MYHCB, CMYP7A, CMYP7B, MYH7
Species	Human
Gene ID (Human)	4625
Gene ID	4625
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1240-1340 of human MYH7 (NP_000248.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P12883
Protein Weight	223 kDa
Shipping	Ice bag