

Anti-Myh7B Rabbit Monoclonal Antibody For Wb, Elisa - A7E2Y1

Item Number: CM0002346

Introduction

Rabbit monoclonal antibody targeting human MYH7B (myosin-7B, MYH14), validated for Western blot and ELISA. Cross-reacts with mouse and rat. Suitable for studies of muscle contraction and cardiac biology.

[Learn More](#)

Parameter	Value
Product Name	MYH7B Rabbit mAb
Remarks/Alias	MHC14; MYH14; IncMYH7b; MYH7B
Species	Human
GeneID (Human)	57644
GeneID	57644
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1750-1983 of human MYH7B (NP_065935.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	A7E2Y1
Protein Weight	226 kDa
Shipping	Ice bag