

Anti-Myh7/ β -Mhc Monoclonal Antibody For Wb, Elisa - P12883

Item Number: CM0006964



Introduction

Rabbit monoclonal antibody targeting human MYH7/ β -MHC (myosin heavy chain 7). Validated for WB and ELISA, with cross-reactivity to mouse and rat. Suitable for studying muscle contraction and cardiomyopathy. Predicted molecular weight 223 kDa.

[Learn More](#)

Parameter	Value
Product Name	MYH7/ β -MHC Rabbit mAb
Remarks/Alias	CMH1; MPD1; SPMD; SPMM; CMD1S; MYHCB; CMYP7A; CMYP7B; MYH7/ β -MHC
Species	Human
Gene ID (Human)	4625
Gene ID	4625
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1700-1800 of human MYH7/ β -MHC (NP_000248.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P12883
Protein Weight	223kDa
Shipping	Ice bag