

Anti-Myh7/ β -Mhc Rabbit Polymab® Monoclonal Antibody For Wb, If-P, Elisa - P12883

Item Number: CM0007099



Introduction

Rabbit recombinant monoclonal antibody against MYH7/ β -MHC (myosin heavy chain 7, cardiac beta isoform). Validated for WB, IF-P, ELISA in human, mouse, rat. Supports research into muscle contraction and cardiomyopathy.

[Learn More](#)

Parameter	Value
Product Name	MYH7/ β -MHC Rabbit PolymAb®
Remarks/Alias	CMH1; MPD1; SPMD; SPMM; CMD1S; MYHCB; CMYP7A; CMYP7B
Species	Human
Gene ID	4625
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1240-1340 of human MYH7/ β -MHC (NP_000248.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P12883
Protein Weight	223kDa
Shipping	Ice bag