


Anti-Myh9 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P35579

Item Number: CM0013875



Introduction

MYH9 Rabbit pAb is a polyclonal antibody validated for WB, IHC-P, and ELISA, reacting with human and mouse MYH9/Myosin-9, crucial for cell motility and cytokinesis research.

[Learn More](#)

Parameter	Value
Product Name	MYH9 Rabbit pAb
Remarks/Alias	MHA; FTNS; EPSTS; BDPLT6; DFNA17; MATINS; NMMHCA; NMHC-II-A; NMMHC-IIA; MYH9
Species	Human
Gene ID (Human)	4627
Gene ID	4627
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1711-1960 of human MYH9 (NP_002464.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	P35579
Protein Weight	227kDa
Shipping	Ice bag