

Anti-Myl1 Polyclonal Antibody For Wb, Elisa - P05976

Item Number: CM0024714

Introduction

High-quality MYL1 rabbit polyclonal antibody validated for WB and ELISA. Reactive with human, mouse, and rat samples. Predicted molecular weight 21 kDa. Ideal for skeletal muscle research.

[Learn More](#)

Parameter	Value
Product Name	MYL1 Rabbit pAb
Remarks/Alias	MLC1; MLC-1; MLC1F; MLC3F; CMYP14; MLC1/3; MYOFTA; MYL1
Species	Human
GeneID (Human)	4632
GeneID	4632
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-150 of human MYL1 (NP_524146.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P05976
Protein Weight	21kDa
Shipping	Ice bag