

Anti-MyI12A Polyclonal Antibody For Wb - P19105

Item Number: CM0025727

Introduction

Rabbit polyclonal antibody targeting human MYL12A (MLC-2B), validated for Western blot and ELISA. Recombinant immunogen spanning amino acids 92-171. Ideal for studying myosin regulation and cell motility.

[Learn More](#)

Parameter	Value
Product Name	MYL12A Rabbit pAb
Remarks/Alias	MLCB, MRCL3, MRLC3, MYL2B, MLC-2B, HEL-S-24, MYL12A
Species	Human
GeneID (Human)	10627
GeneID	10627
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 92-171 of human MYL12A (NP_006462.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P19105
Protein Weight	20kDa
Shipping	Ice bag