

Anti-Myl3 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P08590

Item Number: CM0028959



Introduction

MYL3 Rabbit Polyclonal Antibody (P08590) validated for WB, IF/ICC, ELISA in human, mouse, and rat samples. Ideal for cardiac muscle and myosin research. One-stop sourcing from CamelBio.

[Learn More](#)

Parameter	Value
Product Name	MYL3 Rabbit pAb
Remarks/Alias	CMH8, VLC1, VLCL, MLC1V, MLC1SB, MLC-IV/sb, MYL3
Species	Human
GeneID (Human)	4634
GeneID	4634
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-195 of human MYL3 (NP_000249.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P08590
Protein Weight	22kDa
Shipping	Ice bag