

Anti-Myl6B Rabbit Polyclonal Antibody (Ko Validated) For Wb And Elisa - P14649

Item Number: CM0023137



Introduction

MYL6B (MLC15A) rabbit polyclonal antibody, KO validated, for Western blot and ELISA. Suitable for human, mouse, rat research on myosin light chain function and muscle biology.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] MYL6B Rabbit pAb
Remarks/Alias	MLC15A; 6B
Species	Human
GeneID	140465
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-208 of human MYL6B (NP_002466.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P14649
Protein Weight	23kDa
Shipping	Ice bag