

Anti-Myl9 Monoclonal Antibody For Wb, Elisa - P24844

Item Number: CM0011354

Introduction

Recombinant rabbit monoclonal antibody against MYL9, validated for WB and ELISA. Recognizes human and mouse myosin regulatory light chain 9 (P24844), a key regulator of smooth muscle and nonmuscle contractility. Suitable for cytoskeleton and cell motility studies.

[Learn More](#)

Parameter	Value
Product Name	MYL9 Rabbit mAb
Remarks/Alias	LC20; MLC2; MRLC1; MYRL2; MLC-2C; MMIHS4; MYL9
Species	Human
Gene ID (Human)	10398
Gene ID	10398
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human MYL9 (P24844)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P24844
Protein Weight	20kDa
Shipping	Ice bag