

Anti-Myl9 Rabbit Pab For Wb, If-P, Elisa - P24844

Item Number: CM0030134

Introduction

Rabbit polyclonal antibody against human MYL9 (myosin regulatory light chain 9) for Western blot (WB), immunofluorescence (IF-P), and ELISA. Detects endogenous MYL9 in human, mouse, and rat samples. Useful for smooth muscle and cell motility research.

[Learn More](#)

Parameter	Specification
Product Name	MYL9 Rabbit pAb
Remarks/Alias	LC20; MLC2; MRLC1; MYRL2; MLC-2C; MMIHS4; MYL9
Species	Human
GeneID (human)	10398
GeneID	10398
Immunogen	Synthetic peptide; A synthetic peptide corresponding to a sequence within amino acids 1-100 of human MYL9 (NP_006088.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P24844
Protein Weight	20kDa
Shipping	Ice bag