

Anti-Mylip Polyclonal Antibody For Wb, Elisa - Q8Wy64

Item Number: CM0026708



Introduction

Rabbit polyclonal antibody targeting MYLIP (IDOL/MIR), an E3 ubiquitin ligase critical for cholesterol metabolism. Validated for WB and ELISA, cross-reacts with human, mouse, and rat.

[Learn More](#)

Parameter	Value
Product Name	MYLIP Rabbit pAb
Remarks/Alias	MIR; IDOL; MYLIP
Species	Human
GeneID (human)	29116
GeneID	29116
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 266-445 of human MYLIP (NP_037394.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8WY64
Protein Weight	50kDa
Shipping	Ice bag