

Anti-Myo10 Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Hd67

Item Number: CM0015915



Introduction

Rabbit polyclonal antibody targeting human MYO10 (Myosin-X), suitable for WB, IF/ICC, and ELISA. Involved in filopodia formation and actin-based motility. Cross-reacts with mouse.

[Learn More](#)

Parameter	Value
Product Name	MYO10 Rabbit pAb
Remarks/Alias	MyoX; MYO10
Species	Human
GeneID (Human)	4651
GeneID	4651
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 845-944 of human MYO10 (NP_036466.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9HD67
Protein Weight	237kDa
Shipping	Ice bag