

Anti-Myo18A Rabbit Polyclonal Antibody For Wb And Elisa - Q92614

Item Number: CM0027718

Introduction

Rabbit polyclonal antibody targeting human unconventional myosin-XVIIIa (MYO18A), cross-reactive with mouse and rat. Validated for WB and ELISA, ideal for research in Golgi trafficking, macrophage function, and cytoskeletal dynamics.

[Learn More](#)

Parameter	Value
Product Name	MYO18A Rabbit pAb
Remarks/Alias	MAJN; TIAF1; MYSPDZ; SPR210; SP-R210; MYO18A
Species	Human
Gene ID (Human)	399687
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1970-2054 of human MYO18A (NP_510880.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q92614
Protein Weight	233kDa
Shipping	Ice bag