

Anti-Myo19 Rabbit Monoclonal Antibody For Wb, Elisa - Q96H55

Item Number: CM0003787

Introduction

Rabbit monoclonal antibody against human MYO19 (MYOHD1), validated for WB and ELISA applications. Reacts with human, mouse, rat, and monkey. Suitable for studying mitochondrial dynamics and inheritance.

[Learn More](#)

Parameter	Value
Product Name	MYO19 Rabbit mAb
Remarks/Alias	MYOHD1; MYO19
Species	Human
Gene ID (Human)	80179
Gene ID	80179
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 450-650 of human MYO19 (Q96H55)
Source	Rabbit
Category	Monoclonal Antibody
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat, Monkey
SWISS	Q96H55
Protein Weight	109kDa
Shipping	Ice bag