

Anti-Myo1B Rabbit Monoclonal Antibody For Wb, Elisa - 043795

Item Number: CM0003333

Introduction

Rabbit recombinant monoclonal antibody targeting human MYO1B, validated for WB and ELISA applications. Cross-reacts with human, mouse, and rat samples. Ideal for studying neuronal migration, neurite outgrowth, and vesicular transport.

[Learn More](#)

Parameter	Value
Product Name	MYO1B Rabbit mAb
Remarks/Alias	MMIa, myr1, MYH-1c, MMI-alpha, MYO1B
Species	Human
GeneID (Human)	4430
GeneID	4430
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 916-1136 of Human MYO1B (NP_001123630.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O43795
Protein Weight	132kDa
Shipping	Ice bag