

Anti-Myo1C Polyclonal Antibody For Wb, Ihc-P, Ip, Elisa - 000159

Item Number: CM0024813



Introduction

Rabbit polyclonal antibody against human MYO1C, validated for WB, IHC-P, IP, and ELISA with cross-reactivity to mouse and rat. Targets an unconventional myosin involved in actin-mediated processes and transcription.

[Learn More](#)

Parameter	Value
Product Name	MYO1C Rabbit pAb
Remarks/Alias	NMI; MM1b; myr2; Myo1C; MMI-beta; MYO1C
Species	Human
GeneID (human)	4641
GeneID	4641
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	000159
Protein Weight	122kDa
Shipping	Ice bag