

Anti-Myo1C Rabbit Monoclonal Antibody For Wb, If/Icc, If-P, Ihc-P - O00159

Item Number: CM0007951



Introduction

High-specificity MYO1C Rabbit Monoclonal Antibody validated for WB, IF/ICC, IF-P, IHC-P, ELISA in human, mouse, rat. Ideal for studying actin motor function, insulin signaling, and hair cell adaptation.

[Learn More](#)

Parameter	Value
Product Name	MYO1C Rabbit mAb
Remarks/Alias	NMI; MMIb; myr2; MyoIC; MMI-beta; MYO1C
Species	Human
GeneID (Human)	4641
GeneID	4641
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 769-1028 of human MYO1C (NP_001074248.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O00159
Protein Weight	122kDa
Shipping	Ice bag