


Anti-Myo5B Rabbit Polyclonal Antibody For Wb, If-P, Elisa - Q9ULV0

Item Number: CM0027533



Introduction

Rabbit polyclonal antibody targeting MYO5B, validated for Western blot, immunofluorescence, and ELISA. Reacts with human and mouse. Ideal for vesicular trafficking, cell polarization, and bile acid transport studies.

[Learn More](#)

Parameter	Value
Product Name	MYO5B Rabbit pAb
Remarks/Alias	DIAR2; MVID1; PFIC10; MYO5B
Species	Human
GeneID (human)	4645
GeneID	4645
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1440-1690 of human MYO5B (NP_001073936.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9ULV0
Protein Weight	214kDa
Shipping	Ice bag