

Anti-Myocd Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q8Izq8

Item Number: CM0020051

Introduction

Rabbit polyclonal antibody against human MYOCD (Myocardin), validated for WB, IHC-P, ELISA. Reacts with human and rat. Myocardin is a smooth/cardiac muscle transcriptional cofactor essential for cardiogenesis and myogenesis.

[Learn More](#)

Parameter	Value
Product Name	MYOCD Rabbit pAb
Remarks/Alias	MGBL; MYCD; MYOCD
Species	Human
Gene ID (Human)	93649
Gene ID	93649
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 680-780 of human MYOCD (NP_705832.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Rat
SWISS	Q8IZQ8
Protein Weight	102kDa
Shipping	Ice bag