

Anti-Myof Rabbit Polyclonal Antibody (Ko Validated) For Wb, Elisa - Q9Nzm1

Item Number: CM0018172



Introduction

Rabbit polyclonal antibody against human MYOF (Myoferlin), validated for Western blot and ELISA. Suitable for cross-reactivity with human and mouse samples. Ideal for membrane repair and endocytic recycling studies.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] MYOF Rabbit pAb
Remarks/Alias	HAE7; FER1L3; OF
Species	Human
GeneID (Human)	26509
GeneID	26509
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-190 of human MYOF (NP_038479.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9NZM1
Protein Weight	235kDa
Shipping	Ice bag