


Anti-Myoglobin Monoclonal Antibody For Wb, If-P, Ihc-P, Elisa - P02144

Item Number: CM0006911



Introduction

This rabbit monoclonal antibody against human Myoglobin (MYOSB; PVALB) is validated for WB, IF-P, IHC-P, ELISA and cross-reacts with Mouse, Rat. Ideal for studying oxygen storage, nitrite reductase activity, and muscle tissue biology.

[Learn More](#)

Parameter	Value
Product Name	Myoglobin Rabbit mAb
Remarks/Alias	MYOSB; PVALB; Myoglobin
Species	Human
Gene ID (Human)	4151
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 2-154 of human Myoglobin(NP_005359.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P02144
Protein Weight	17kDa
Shipping	Ice bag