


Anti-Myoglobin Rabbit Polyclonal Antibody For Wb, If-P, Elisa - P02144

Item Number: CM0026272



Introduction

High-quality rabbit polyclonal antibody targeting human myoglobin, validated for Western blot, IF-P, and ELISA. Cross-reacts with mouse and rat. Ideal for oxygen storage and muscle tissue research.

[Learn More](#)

Parameter	Value
Product Name	Myoglobin Rabbit pAb
Remarks/Alias	MYOSB; PVALB; Myoglobin
Species	Human
Gene ID (Human)	4151
Gene ID	4151
Immunogen	Recombinant protein - Recombinant fusion protein containing a sequence corresponding to amino acids 1-154 of human Myoglobin (NP_005359.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P02144
Protein Weight	17kDa
Shipping	Ice bag