

Anti-Myosin Heavy Chain Rabbit Monoclonal Antibody For Wb, Ihc, If, Elisa - P12883 / P13533

Item Number: CM0009434



Introduction

Rabbit monoclonal antibody targeting human Myosin heavy chain (MYH7/MYH6), validated in WB, IHC-P, IF-P, ELISA. Reacts with mouse, rat. Ideal for muscle research and IVD assays.

[Learn More](#)

Parameter	Specification
Product Name	Myosin heavy chain Rabbit mAb
Remarks/Alias	CMD1S; CMH1; MPD1; MYHCB; SPMD; SPMM; Myosin heavy chain
Species	Human
GeneID (Human)	4625/4624
GeneID	4625/4624
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P12883/P13533
Protein Weight	224 kDa
Shipping	Ice bag