

Anti-Myotubularin (Mtm1) Polyclonal Antibody For Wb, Elisa - Q13496

Item Number: CM0017088

Introduction

Rabbit polyclonal antibody against human myotubularin (MTM1) for WB and ELISA. Reacts with human, mouse, rat. Useful for studying phosphoinositide metabolism, endosomal trafficking, and muscle diseases.

[Learn More](#)

Parameter	Value
Product Name	MTM1 Rabbit pAb
Remarks/Alias	CNM; CNMX; MTMX; XLMTM; MTM1
Species	Human
Gene ID (Human)	4534
Gene ID	4534
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 484-603 of human MTM1 (NP_000243.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13496
Protein Weight	70kDa
Shipping	Ice bag