

Anti-Mypn Rabbit Polyclonal Antibody For Wb, Elisa - Q86Tc9

Item Number: CM0022286

Introduction

This rabbit polyclonal antibody against MYPN (myopalladin) is validated for Western blot (WB) and ELISA, detecting human and mouse MYPN. Suitable for sarcomere structure and muscle disease research.

[Learn More](#)

Parameter	Value
Product Name	MYPN Rabbit pAb
Remarks/Alias	MYOP; RCM4; CMH22; NEM11; CMD1DD; MYPN
Species	Human
GeneID (human)	84665
GeneID	84665
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-250 of human MYPN (NP_115967.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q86TC9
Protein Weight	145kDa
Shipping	Ice bag