


Anti-Myrf Polyclonal Antibody For Wb And Elisa - Q9Y2G1

Item Number: CM0023853



Introduction

Anti-MYRF rabbit polyclonal antibody validated for WB and ELISA. Reacts with human and mouse MYRF, a key transcription factor in CNS myelination and oligodendrocyte maturation.

[Learn More](#)

Parameter	Value
Product Name	MYRF Rabbit pAb
Remarks/Alias	MRF; CUGS; MMERV; Ndt80; 11orf9; pqn-47; C11orf9; MYRF
Species	Human
GeneID (human)	745
GeneID	745
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human MYRF (NP_001120864.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9Y2G1
Protein Weight	124kDa
Shipping	Ice bag