

Anti-Mysm1 Rabbit Polyclonal Antibody For Wb, Elisa - Q5Vvj2

Item Number: CM0024810



Introduction

MYSM1 Rabbit Polyclonal Antibody validated for WB and ELISA. Reacts with human, mouse, and rat. Ideal for deubiquitinase research, hematopoietic stem cell function, and immune response studies.

[Learn More](#)

Parameter	Value
Product Name	MYSM1 Rabbit pAb
Remarks/Alias	2ADUB; BMFS4; 2A-DUB; MYSM1
Species	Human
GeneID (Human)	114803
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 720-821 of human MYSM1 (NP_001078956.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q5VVJ2
Protein Weight	95kDa
Shipping	Ice bag