


Anti-Myt1 Rabbit Polyclonal Antibody For Wb, Elisa - Q01538

Item Number: CM0012971



Introduction

High-quality anti-MYT1 rabbit polyclonal antibody validated for WB and ELISA. Reacts with human, mouse, rat MYT1, a transcription factor crucial for oligodendrocyte development. Ideal for neuroscience and developmental studies.

[Learn More](#)

Parameter	Value
Product Name	MYT1 Rabbit Polyclonal Antibody
Remarks/Alias	MTF1; MYTI; NZF2; PLPB1; ZC2H2C1; ZC2HC4A; C20orf36; MYT1
Species	Human
GeneID (Human)	4661
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 650-820 of human MYT1 (NP_004526.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q01538
Protein Weight	122kDa
Shipping	Ice bag