


Anti-Mzb1 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q8Wu39

Item Number: CM0014022



Introduction

High-quality rabbit polyclonal antibody targeting human MZB1 (pERp1, MEDA-7) for Western blot, IHC-P, and ELISA. Cross-reacts with mouse and rat. Ideal for studying ER protein chaperone functions and B-cell biology.

[Learn More](#)

Parameter	Value
Product Name	MZB1 Rabbit pAb
Remarks/Alias	PACAP; pERp1; MEDA-7; MZB1
Species	Human
GeneID (Human)	51237
GeneID	51237
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-189 of human MZB1 (NP_057543.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8WU39
Protein Weight	21kDa
Shipping	Ice bag