

Anti-Mzf1 Polyclonal Antibody For Wb, Elisa - P28698

Item Number: CM0017496



Introduction

MZF1 Rabbit Polyclonal Antibody for WB and ELISA, validated in human, mouse, rat. Recombinant protein immunogen targeting myeloid zinc finger 1 transcription factor involved in hematopoietic differentiation.

[Learn More](#)

Parameter	Value
Product Name	MZF1 Rabbit pAb
Remarks/Alias	MZF-1; MZF1B; ZFP98; ZNF42; ZSCAN6; MZF1
Species	Human
Gene ID (Human)	7593
Immunogen	Recombinant protein: Recombinant fusion protein containing a sequence corresponding to amino acids 100-350 of human MZF1 (NP_003413.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P28698
Protein Weight	82kDa
Shipping	Ice bag