


Anti-Mzt2B Polyclonal Antibody For Wb, Elisa - Q6Nz67

Item Number: CM0027445



Introduction

Rabbit polyclonal antibody targeting human MZT2B (FAM128B), validated for WB and ELISA. Ideal for studying centrosome function, spindle assembly, and gamma-tubulin ring complex dynamics.

[Learn More](#)

Parameter	Value
Product Name	MZT2B Rabbit pAb
Remarks/Alias	FAM128B; MOZART2B; MZT2B
Species	Human
GeneID (Human)	80097
GeneID	80097
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-158 of human MZT2B (NP_079305.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q6NZ67
Protein Weight	16kDa
Shipping	Ice bag