

Anti-N-Cadherin Polyclonal Antibody For Wb, If-P, Ihc-P, Elisa - P19022

Item Number: CM0029385

Introduction

Rabbit polyclonal antibody against human N-Cadherin, validated for WB, IF-P, IHC-P, ELISA. Cross-reacts with mouse and rat. Suitable for cell adhesion and neuroscience studies.

[Learn More](#)

Parameter	Value
Product Name	N-Cadherin Rabbit pAb
Remarks/Alias	CDHN; NCAD; ACOGS; ADHD8; CD325; ARVD14; CDw325; N-Cadherin
Species	Human
GeneID (Human)	1000
GeneID	1000
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 488-587 of human N-Cadherin (NP_001783.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P19022
Protein Weight	100kDa
Shipping	Ice bag