

Anti-N-Cadherin Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P19022

Item Number: CM0014183

Introduction

Rabbit polyclonal antibody against human N-Cadherin (CDH2, Neural cadherin), validated for WB, IHC-P, ELISA. Cross-reacts with mouse and rat. Suitable for cell adhesion, neurobiology, and cancer research.

[Learn More](#)

Parameter	Value
Product Name	N-Cadherin Rabbit pAb
Remarks/Alias	CDHN; NCAD; ACOGS; ADHD8; CD325; ARVD14; CDw325; N-Cadherin
Species	Human
Gene ID (Human)	1000
Gene ID	1000
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 807-906 of human N-Cadherin (NP_001783.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB,IHC-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P19022
Protein Weight	100kDa
Shipping	Ice bag