


Anti-N-Wasp/Wasl Monoclonal Antibody For Wb, Elisa - 000401

Item Number: CM0006805



Introduction

Validated rabbit monoclonal antibody specific to human N-WASP/WASL (000401). Detects a 55kDa band in WB and is suitable for ELISA. Predicted cross-reactivity with mouse and rat. Ideal for actin cytoskeleton and cell motility studies.

[Learn More](#)

Parameter	Value
Product Name	N-WASP/WASL Rabbit mAb
Remarks/Alias	NWASP; WASPB; N-WASP; N-WASP/WASL
Species	Human
GeneID (Human)	8976
GeneID	8976
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 350-450 of human N-WASP/WASL (000401).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O00401
Protein Weight	55kDa
Shipping	Ice bag