


Anti-N-Wasp/Wasl Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - O00401

Item Number: CM0028159



Introduction

Rabbit polyclonal antibody targeting human N-WASP/WASL. Validated for WB, IHC-P, IF/ICC, and ELISA. Cross-reacts with human, mouse, and rat. Recombinant immunogen (aa 1-250). Ideal for actin nucleation and cytoskeleton studies.

[Learn More](#)

Parameter	Value
Product Name	N-WASP/WASL Rabbit pAb
Remarks/Alias	NWASP; WASPB; N-WASP; N-WASP/WASL
Species	Human
GeneID (Human)	8976
GeneID	8976
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-250 of human N-WASP/WASL (NP_003932.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O00401
Protein Weight	55kDa
Shipping	Ice bag