


Anti-N4BP1 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - O75113

Item Number: CM0026329



Introduction

High-quality rabbit polyclonal antibody targeting N4BP1 (NEDD4-binding protein 1), validated for WB, IF/ICC, and ELISA. Recognizes human and mouse N4BP1. Ideal for innate immunity and viral restriction research.

[Learn More](#)

Parameter	Value
Product Name	N4BP1 Rabbit pAb
Remarks/Alias	N4BP1
Species	Mouse
GeneID (human)	9683
GeneID	80750
Immunogen	Recombinant protein; Recombinant fusion protein containing a sequence corresponding to amino acids 190-590 of mouse N4BP1 (NP_085040.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	O75113
Protein Weight	100kDa
Shipping	Ice bag