

Anti-N6Amt1/Hemk2 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Y5N5

Item Number: CM0028414



Introduction

Rabbit polyclonal antibody targeting human N6AMT1/HEMK2 (KMT9) validated for WB, IF/ICC, and ELISA; cross-reacts with mouse and rat. Ideal for chromatin regulation and methyltransferase studies.

[Learn More](#)

Parameter	Value
Product Name	N6AMT1/HEMK2 Rabbit pAb
Remarks/Alias	KMT9; MTQ2; PrmC; HEMK2; N6AMT; PRED28; C21orf127; m.HsaHemK2P; N6AMT1/HEMK2
Species	Human
GeneID (Human)	29104
GeneID	29104
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-186 of human N6AMT1/HEMK2 (NP_877426.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y5N5
Protein Weight	23kDa
Shipping	Ice bag