

Anti-N6Amt1/Hemk2 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9Y5N5

Item Number: CM0018983



Introduction

Rabbit polyclonal antibody against human N6AMT1/HEMK2 methyltransferase, validated for WB, IHC-P, ELISA. Cross-reacts with human and mouse; suitable for chromatin regulation and arsenic toxicity research.

[Learn More](#)

Parameter	Value
Product Name	N6AMT1/HEMK2 Rabbit pAb
Remarks/Alias	KMT9; MTQ2; PrmC; HEMK2; N6AMT; PRED28; C21orf127; m.HsaHemK2P; N6AMT1/HEMK2
Species	Human
Gene ID	29104
Immunogen	Recombinant protein
Source	rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9Y5N5
Protein Weight	23kDa
Shipping	Ice bag