

Anti-Naa10 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P41227

Item Number: CM0028477



Introduction

High-quality anti-NAA10 rabbit polyclonal antibody validated for WB, IHC-P, IF/ICC, and ELISA. Targets human NAA10 with cross-reactivity to mouse and rat. Ideal for acetylation research and cancer studies.

[Learn More](#)

Parameter	Value
Product Name	NAA10 Rabbit pAb
Remarks/Alias	MAA; TE2; ARD1; LZMS; NATD; ARD1A; ARD1P; OGDNS; hARD1; DXS707; MCOP51; NAA10
Species	Human
GeneID (Human)	8260
GeneID	8260
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 161-235 of human NAA10 (NP_003482.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P41227
Protein Weight	26kDa
Shipping	Ice bag