


Anti-Naa20 Rabbit Polyclonal Antibody For Wb, Ihc-P, Ip, Elisa - P61599

Item Number: CM0030811



Introduction

NAA20 Rabbit pAb is a polyclonal antibody for detection of N-alpha-acetyltransferase 20. Suitable for WB, IHC-P, IP, ELISA with cross-reactivity to human, mouse, and rat. Useful for N-terminal acetylation, cell migration, and actomyosin fiber studies.

[Learn More](#)

Parameter	Value
Product Name	NAA20 Rabbit pAb
Remarks/Alias	NAT3; NAT5; MRT73; NAT3P; NAT5P; dj1002M8.1; NAA20
Species	Human
Gene ID (Human)	51126
Gene ID	51126
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 49-178 of human NAA20 (NP_057184.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P61599
Protein Weight	20kDa
Shipping	Ice bag