

Anti-Naa50 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Gzz1

Item Number: CM0026289



Introduction

N-alpha-acetyltransferase 50 (NAA50) Rabbit Polyclonal Antibody validated for WB, IF/ICC, ELISA. Cross-reacts with human, mouse, rat. Suitable for studies of protein N-terminal acetylation, chromatin cohesion, and epigenetic regulation.

[Learn More](#)

Parameter	Value
Product Name	NAA50 Rabbit pAb
Remarks/Alias	SAN; MAK3; NAT5; NAT13; NAT5P; NAT13P; hNaa50p; NAA50
Species	Human
Gene ID	80218
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-169 of human NAA50 (NP_079422.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9GZZ1
Protein Weight	19kDa
Shipping	Ice bag