

Anti-Naa50 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q9Gzz1

Item Number: CM0028401



Introduction

Rabbit polyclonal antibody targeting human NAA50 (Q9GZZ1), validated for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with human, mouse, and rat samples. Ideal for studying N-terminal acetylation and sister chromatid cohesion.

[Learn More](#)

Parameter	Value
Product Name	NAA50 Rabbit pAb
Remarks/Alias	SAN; MAK3; NAT5; NAT13; NAT5P; NAT13P; hNaa50p; NAA50
Species	Human
Gene ID	80218
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-169 of human NAA50 (NP_079422.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9GZZ1
Protein Weight	19kDa
Shipping	Ice bag