


## Anti-Nab1 Polyclonal Antibody For Wb, Elisa - Q13506

Item Number: CM0019102



### Introduction

High-quality rabbit polyclonal antibody against human NAB1 (NGFI-A-binding protein 1), validated for WB and ELISA. Detects endogenous NAB1, a transcriptional repressor of EGR1 and EGR2. Cross-reacts with mouse.

[Learn More](#)

| Parameter        | Value                                                                                                             |
|------------------|-------------------------------------------------------------------------------------------------------------------|
| Product Name     | NAB1 Rabbit pAb                                                                                                   |
| Remarks/Alias    | NAB1                                                                                                              |
| Species          | Human                                                                                                             |
| GeneID (Human)   | 4664                                                                                                              |
| GeneID           | 4664                                                                                                              |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 189-317 of human NAB1 (NP_005957.2) |
| Source           | Rabbit                                                                                                            |
| Category         | Polyclonal Antibodies                                                                                             |
| Application      | WB, ELISA                                                                                                         |
| Cross Reactivity | Human, Mouse                                                                                                      |
| SWISS            | Q13506                                                                                                            |
| Protein Weight   | 54kDa                                                                                                             |
| Shipping         | Ice bag                                                                                                           |