

Anti-Nab2 Rabbit Polyclonal Antibody For Wb, Elisa - Q15742

Item Number: CM0027528

Introduction

High-quality anti-NAB2 rabbit polyclonal antibody validated for WB and ELISA. Targets human and mouse NGFI-A-binding protein 2, a transcriptional repressor of EGR1/EGR2. Ideal for transcription regulation studies.

[Learn More](#)

Parameter	Value
Product Name	NAB2 Rabbit pAb
Remarks/Alias	MADER; NAB2
Species	Human
Gene ID (Human)	4665
Gene ID	4665
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human NAB2 (NP_005958.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q15742
Protein Weight	57kDa
Shipping	Ice bag