


Anti-Naca Polyclonal Antibody For Wb, Elisa - Q13765

Item Number: CM0015695



Introduction

Rabbit polyclonal antibody against human NACA (NAC-alpha) for WB and ELISA. Cross-reacts with mouse and rat. Target is nascent polypeptide-associated complex subunit alpha involved in protein targeting and transcription.

[Learn More](#)

Parameter	Value
Product Name	NACA Rabbit pAb
Remarks/Alias	HSD48; NACA1; skNAC; NAC-alpha; NACA
Species	Human
Gene ID (Human)	4666
Gene ID	4666
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13765
Protein Weight	23kDa
Shipping	Ice bag