


Anti-Naca Rabbit Monoclonal Antibody For Wb, Elisa - Q13765

Item Number: CM0007238



Introduction

Recombinant rabbit monoclonal antibody targeting human NACA (Alpha-NAC, nascent polypeptide-associated complex subunit alpha). Validated for WB and ELISA applications; cross-reacts with mouse and rat. Useful for ER targeting and gene regulation studies.

[Learn More](#)

Parameter	Value
Product Name	NACA Rabbit mAb
Remarks/Alias	HSD48; NACA1; skNAC; NAC-alpha; NACA
Species	Human
GeneID (Human)	4666
GeneID	4666
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-215 of human NACA (NP_005585.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13765
Protein Weight	23kDa
Shipping	Ice bag